

GUIDE TO SERVICES

TREATMENTS

- Hemofiltration On-line pre dilution, post/mixed
- Bicarbonate dialysis with very biocompatible membranes
- Single needle dialysis

MONITORING

- Sodium profiles and ultrafiltration
- OCM, dialytic adequacy constant monitoring
- BTM, temperature and recirculation management form
- BPM, arterial pressure constant monitoring
- BCM, corporal composition monitoring and nutritional status
- BVM, Monitoring of blood volume changes and Ultrafiltration

EMODIALYSIS MONITORS

- Cordiax 5008

DYALISIS PLACES

- n°13 technical workstations for hemodialysis

	Monday	Tuesday	Wednesday	Thursday	Friday	Saturday	Sunday
OPENING DAYS AND HOURS	07.00-23.00	07.00-19.00	07.00-23.00	07.00-19.00	07.00-23.00	07.00-19.00	<i>closed</i>
DIALYSIS SHIFTS	n°3 (M-A-E)	n°2 (M-A)	n°3 (M-A-E)	n°2 (M-A)	n°3 (M-A-E)	n°2 (M-A)	-

M = morning | A = afternoon | E = evening | For timetable details, please contact the dialysis center

STAFF AND SERVICES

- Nephrologist Responsible Director: **Dott. Raffaele Galato**
E-mail: Raffaele.Galato@freseniusmedicalcare.com
- Dialysis session conductor doctors
- Professional nurses
- S.S.N. Accreditate center**
- Charges to be borne by S.S.N.**

